

### Compression Levels

### Clinical Indications

**Circular Knit**  
12-16 mmHg  
15-20 mmHg

- aching/fatigued legs
- prophylaxis during pregnancy
- mild edema of ankle and foot
- mild varicosities
- post sclerotherapy

**Circular Knit**  
20-30 mmHg  
**Flat Knit**  
18-21 mmHg

- aching/fatigued legs
- prophylaxis during pregnancy
- hereditary tendency toward varicose veins
- mild lymphedema
- hypertrophic scar treatment
- prevention/recurrence of venous ulcers
- mild venous insufficiency
- moderate varicosities
- post sclerotherapy
- burn scar management
- DVT prevention

**Circular Knit**  
30-40 mmHg  
**Flat Knit**  
23-32 mmHg

- chronic venous insufficiency
- severe varicosities
- post surgical
- post-phlebectomy
- burn scar management
- moderate lymphedema
- prevention/recurrence of venous ulcers
- management of active venous ulcers
- moderate and post traumatic edema
- post sclerotherapy
- pronounced varicosities during pregnancy
- orthostatic hypotension
- prevention of DVT/post thrombotic syndrome

**Circular Knit**  
40-50 mmHg  
**Flat Knit**  
34-46 mmHg

- management of active venous ulcers
- severe chronic venous insufficiency
- severe tendency toward edema
- severe lymphedema

**Flat Knit**  
50+ mmHg

- severe post-thrombotic conditions
- severe lymphedema
- elephantiasis

### Contra indications

- open or untreated ulcers
- acute thrombophlebitis
- arterial disease
- uncontrolled congestive heart failure
- weeping dermatosis
- phlebothrombosis
- acute dermatitis
- intermittent claudication










For more information or to locate a Juzo dealer, contact us at:

Phone: 1 (800) 222-4999

Fax: 1 (800) 645-2519

[www.juzousa.com](http://www.juzousa.com)

[support@juzousa.com](mailto:support@juzousa.com)

Product	Extremity	Stages	Model Number
<b>Juzo Soft 2000</b> <ul style="list-style-type: none"> <li>extremely soft</li> <li>inconspicuous when worn</li> </ul>		<b>Latency Stage 1</b>	<b>2001 (20-30 mmHg) 2002 (30-40 mmHg)</b>
<b>Juzo Varin Soft-in</b> <ul style="list-style-type: none"> <li>soft inside for greater ease in application</li> <li>ideal for sensitive skin</li> <li>durable, opaque knit</li> </ul>		<b>Stage 1 stage 2</b>	<b>3511 (20-30 mmHg) 3512 (30-40 mmHg) 3513 (40-50 mmHg)</b>
<b>Juzo Silver Varin Soft-in</b> <ul style="list-style-type: none"> <li>antimicrobial, anti-static thermodynamic</li> <li>ideal for patients with skin breakdown</li> </ul>		<b>Stage 1 Stage 2</b>	<b>3561 (20-30 mmHg) 3562 (30-40mmHg) 3563 (40-50 mmHg)</b>
<b>Juzo Varin</b> <ul style="list-style-type: none"> <li>durable, opaque knit</li> <li>wide assortment of styles and sizes</li> </ul>		<b>Stage 1 Stage 2</b>	<b>3511 (20-30 mmHg) 3512 (30-40 mmHg) 3513 (40-50 mmHg)</b>
<b>Juzo Varin Full Knit</b> <ul style="list-style-type: none"> <li>durable, opaque knit</li> <li>continued compression in thighs and abdomen</li> </ul>		<b>Stage 2</b>	<b>3512L (30-40 mmHg) 3513L (40-50 mmHg)</b>
<b>Juzo Helastic</b> <ul style="list-style-type: none"> <li>dense knit for optimal therapeutic effect</li> <li>can be made in any shape or size</li> </ul>		<b>Stage 2 Stage 3</b>	<b>3021 (18-21 mmHg) 3022 (23-32 mmHg) 3023 (34-46 mmHg) 3024 (50+ mmHg)</b>
<b>Juzo Strong</b> <ul style="list-style-type: none"> <li>micro-massaging weave</li> <li>maximum containment</li> <li>can be made in any shape or size</li> </ul>		<b>Stage 2 Stage 3</b>	<b>3051 (18-21 mmHg) 3052 (23-32 mmHg) 3053 (34-46 mmHg) 3054 (50+mmHg)</b>
<b>Juzo Silver Helastic</b> <ul style="list-style-type: none"> <li>Flat knit made with silver</li> <li>Anti-odor, antimicrobial, anti static</li> <li>Machine wash and dry</li> </ul>		<b>Stage 2 Stage 3</b>	<b>3021SV (18-21 mmHg) 3022SV (23-32 mmHg) 3023SV (34-46 mmHg) 3024SV (50+mmHg)</b>
<b>Juzo Silver Strong</b> <ul style="list-style-type: none"> <li>micro-massaging weave</li> <li>maximum containment</li> <li>Eliminates 99% of Staph, Strep and Candida Albicans</li> </ul>		<b>Stage 2 Stage 3</b>	<b>3051SV (18-21 mmHg) 3052SV (23-32 mmHg) 3053SV (34-46 mmHg) 3054SV (50+mmHg)</b>